Standard Operating Procedures

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Title	세포계수기		
SOP No.	IAI/EQM/032	Total page	11
Version	01	Effective date	

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Instrument Type	Bench-top cell counter
Dimensions (WxDxH)	22 cm x 21 cm X 9 cm
Weight	1.8 kg (without the power cord/adapter)
Light Source	LED
Excitation Wave Length	470 +/- 20 nm (Blue)
Emission Wave Length	525 +/- 25 nm (Green), 600 LP (Red)
Cell Concentration Range	$1 \times 10^4 \sim 1 \times 10^7$ (optimally $5 \times 10^4 \sim 1 \times 10^7$) cells/ml
Cell Diameter Range	1~90 (optimally 5~60) μm
Cell Circularity Range	30~100%
Cell Viability Range	0~100%
Image Resolution	CMOS camera (5 MP)
Image Type	TIF
Reporting	PDF report, .CSV file
LCD Display	7 inches (800 x 480 pixels)
Processing Time	Approx. 7 seconds for bright field & approx. 30 seconds for dual fluorescence (may depend on cell type and concentration)

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
(3) Report : PDF report 저장여부

4.2.2.2. Counting Options


(1) With Trypan Blue : Trypan blue 염색 이용시 on

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Dilution Calculator

Current Concentration
4.5
× 10^e
5
/mL
of
Total

Desired Concentration
× 10^e
/mL
Final Volume
mL

1 2 3 4 5 6 7 8 9 0

. Calculate

☒ Analyzed Image
☒ Raw Image
☒ Report
Print

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
Problem	Possible Cause	Solution
Inaccurate cell count	Unfocused cell image	Ensure that the cell image is appropriately focus by using the focus knob. You may use the "Zoom-in" button when adjusting the cell image. Make sure that live cells have bright centers and dead cells have dark/blue centers.
	Clumped cells	Make sure that cells are not clumped. The more single cells, the better counting results.
	Cell concentration range	The cell concentration range of the instrument is preferably $5 \times 10^4 \sim 1 \times 10^7$ cells/ml. Make sure that the sample is in this range. The sample may need to be concentrated or diluted.
	Counting slide insertion	Ensure that the counting slide is inserted completely into the instrument. When the slide is inserted, a soft click can be heard.
	Sample loading	If the counting slide is over- or under-loaded with the sample, it may affect counting results. The optimal amount of sample is $10 \sim 12 \mu\text{l}$ of cell suspension.
	Malfunction of optical components	Any of the optical components may be damaged. Or, the objective lens may be dirty due to dust, spilled samples, or unknown causes. Please contact your local supplier.
	Damaged counting surface	Make sure that the counting area of the slide is transparent before loading the sample. Wear gloves while handling the slide.
	Floating cells	After loading the sample in the counting chip, wait for about 1 min to allow the cells to settle down. Time required for each cell type can vary and should be determined empirically.
	Omitted dyes	In the "Fluorescence Cell Counting" mode, cells cannot be counted without fluorescent nucleic acid binding dyes, e.g., AO and PI.
	Too low or too high fluorescence setting	Make sure fluorescence settings in the preview screen are in a reasonable range.

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